

## Our line of durable Sampling Scoops designed for easy sample retrieval of granules, powders and pastes

Available in multiple sizes and intended for single use to prevent contamination of bulk sample material

### Durable Design

- Large, wide openings for easy sample retrieval
- Flat bottom to stabilize sample on scale or flat surface
- Lot traceable
- Shatterproof constructed of inert, BSE/TSE free material

### Safe Sample Recovery

- Disposable to ensure no chemical cross-contamination
- Made from FDA food grade polystyrene in an ISO 8 Class clean room
- Meets quality control and GMP requirements for sampling
- Available in both, sterile\* and non-sterile for intended usage


\*Electron beam sterilized and individually wrapped

Certificate of Irradiation available upon request



Manufactured  
and packaged  
in a clean  
room

STERILE R



60 mL

125 mL

250 mL

500 mL



### Small Sampling Scoop

Item No.	Description	Volume	L x W x H <i>in</i>	L x W x H <i>mm</i>	Sterile	UOM
120858	Small Sampling Scoop Non-Sterile	60 mL/2 oz.	8.1 x 1.6 x 1.1	205 x 40 x 27		Case 200; 20 x 10/Bulk
120857	Small Sampling Scoop Sterile	60 mL/2 oz.	8.1 x 1.6 x 1.1	205 x 40 x 27	✓	Case 200; 20 x 10/Single Bag



### Medium Sampling Scoop

Item No.	Description	Volume	L x W x H <i>in</i>	L x W x H <i>mm</i>	Sterile	UOM
120856	Medium Sampling Scoop Non-Sterile	125 mL/4 oz.	8.6 x 1.9 x 1.3	218 x 49 x 38		Case 200; 20 x 10/Bulk
120855	Medium Sampling Scoop Sterile	125 mL/4 oz.	8.6 x 1.9 x 1.3	218 x 49 x 38	✓	Case 200; 20 x 10/Single Bag



### Large Sampling Scoop

Item No.	Description	Volume	L x W x H <i>in</i>	L x W x H <i>mm</i>	Sterile	UOM
120854	Large Sampling Scoop Non-Sterile	250 mL/8 oz.	7.7 x 2.3 x 1.7	195 x 58 x 44		Case 200; 20 x 10/Bulk
120853	Large Sampling Scoop Sterile	250 mL/8 oz.	7.7 x 2.3 x 1.7	195 x 58 x 44	✓	Case 200; 20 x 10/Single Bag



### XLarge Sampling Scoop

Item No.	Description	Volume	L x W x H <i>in</i>	L x W x H <i>mm</i>	Sterile	UOM
120852	XLarge Sampling Scoop Non-Sterile	500 mL/16 oz.	10.0 x 3.6 x 2.1	255 x 92 x 53		Case 200; 20 x 10/Bulk
120851	XLarge Sampling Scoop Sterile	500 mL/16 oz.	10.0 x 3.6 x 2.1	255 x 92 x 53	✓	Case 200; 20 x 10/Single Bag