




Aluminum Weighing Boats

Boats for general-purpose weighing, dispensing, storage, food processing, evaporating or drying

Aluminum Weighing Boats

- Sturdy yet lightweight boats feature a smooth, flat bottom for sample handling
- Finger tab handle can be etched for identification prior to weighing





Aluminum Weighing Boats

Item No.	Volume <i>mL</i>	Dimensions <i>in</i>	Dimensions <i>mm</i>	Packaging	Tab	Wall Finish
120937	140 mL	4.1 x 0.8	103 x 20	100		Smooth
120936	80 mL	2.9 x 0.7	73 x 18	100		Smooth
120930	70 mL	2.8 x 0.7	71 x 17	100		Crinkle
HS14521C	50 mL	2.8 x 0.6	70 x 16	100	x	Crinkle
120935	42 mL	2.2 x 0.6	57 x 16	100		Smooth
HS14521B	34 mL	2.2 x 0.6	57 x 14	100	x	Crinkle
120931	25 mL	2 x 0.5	52 x 13	100		Crinkle
HS14521A	20 mL	1.7 x 0.5	43 x 12	100	x	Crinkle
120932	5 mL	1.1 x 0.4	28 x 10	100	x	Crinkle