


Dishes for general-purpose weighing, dispensing, storage, food processing, evaporating or drying

Aluminum Weighing Dish

- Designed for moisture determination balances with slightly tapered walls for stacking

Aluminum Weighing Dish



120933
1 mL



120940
0.25 mL

Aluminum Weighing Dish

Item No.	Volume mL	Dimensions in	Dimensions mm	Packaging	Tab
120938	45 mL	3.5 x 0.3	90 x 7	50	
HS14522	40 mL	4 x 0.3	101 x 8	50	
120934	15 mL	2.7 x 0.2	70 x 6	50	
120933	1 mL	0.98 x 0.3	25 x 7	50	
120940	0.25 mL	0.5 x 0.1	13 x 3.5	100	x