

## The Perfect Grip for Optimal Precision

These two styles of hand-operated grippers are designed to pick up cryogenic vials and crimped vaccine/serum vials by the cap, respectively.

Made to grab with precision and safely remove tubes and vials from cold temperature storage boxes.

Both grippers are crafted from autoclavable, reinforced polypropylene, and feature an ergonomic squeeze area for optimal precision.

Each style is available in a pack of 5 assorted colors for easy identification and organization. A combo pack with one of each style is also available.



Vaccine Vial Gripper head designed fit for vials up to 21 mm diameter

Marked thumb area for optimal ergonomics and control

Cryo Vial Gripper head designed fit for vials up to 13.5 mm diameter



Improved Dexterity while handling vials with thick cryogloves



Slender Cryo Gripper Design handles vials in tight spaces with ease





Ridged Vaccine Gripper Head creates the perfect grip on crimped vaccine/serum vials

Available in 5 colors for easy identification & organization




PK 5 / CS 100

Item No.	Gripper Type	Color	L x W x H in	L x W x H cm
120908	Cryogenic	Assorted* 	5.6 x 0.6 x 0.9	14.3 x 1.5 x 2.3
120909	Vaccine	Assorted* 	6 x 0.8 x 1.4	15.2 x 2 x 3.5

Assorted colors: red, blue, green, purple, yellow.

PK 2 / CS 200

Item No.	Gripper Type	Color	L x W x H in	L x W x H cm
120910	Combination	Blue 	6 x 0.8 x 1.4 (Vaccine) 5.6 x 0.6 x 0.9 (Cryo)	15.2 x 2 x 3.5 (Vaccine) 14.3 x 1.5 x 2.3 (Cryo)

