

Specimen or Sample Transport Box

Ideal for transporting biological and clinical samples, as well as instruments and products that you want to keep clean and dry under testing environmental conditions

The DuraPorter® transparent design is OSHA 29 CFR 1910.1031 and UN3373 compliant

- Polycarbonate construction and available in three colors for user/owner identification
- Accommodates 13 mm, 16 mm and 20 mm tubes with our Half and Full Size OneRack®
- Use the separators provided to create your own compartment sizes



| Item No. | Color | | L x W x H in | L x W x H cm |
|----------|---------------------------|---|-----------------|-----------------|
| 120052 | Clear with Blue Handle | ■ | | |
| 120077 | Red with Red Handle | ■ | 15 x 7.7 x 6.3 | 38 x 19.6 x 16 |
| 120078 | Yellow with Yellow Handle | ■ | | |





Load



Lock



Transport



Store
can be placed in a freezer or refrigerator to keep samples chilled

Lid Opens a Full 180°
to expedite loading and emptying



Extra Security Feature
Accommodates lockout tags (not included)



Included separators can also be used to create inner compartments.
Or, use without racks or separators as one large storage space for larger samples or supplies

A Perfect Fit
OneRack® Tube Rack
Accommodates 13 mm, 16 mm and 20 mm tubes with our Half and Full Size OneRack®

