


Economical, Reusable & Moisture Resistant Plastic

The freezer boxes by Heathrow Scientific, were designed with sustainability in mind and are made from a specially formulated corrugated or thin film polypropylene


With superior durability, polypropylene outlasts the cardboard box alternatives that are susceptible to softening and producing dust and can freeze (down to temperatures of -80°C) and thaw multiple times without an issue

Easy to assemble boxes made from polypropylene material

- Moisture resistant, more durable than cardboard boxes
- Dividers included with the boxes
- Corrugated Polypropylene: Store down to temperatures of -80°C with appropriate temperature viable tubes
- Thin Film: Can be used in mechanical freezers & vessels down to -196°C with appropriate temperature viable tubes
- Vapor phase only



Freezer boxes are supplied "flat-packed"
Eliminating the hassle of stacking and reduces the amount of space needed for storage





Corrugated Polypropylene 0.2/0.5/1.5–2.0/5 mL/Vial



Pack of 10

Item No.	Color		Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120379	White	<input type="checkbox"/>	0.2	144	7.9	5.2 x 5.1 x 1.2	13.3 x 13 x 3
120383	White	<input type="checkbox"/>	0.5	81	12	5.2 x 5.1 x 1.6	13.3 x 13 x 4.2
120262	White	<input type="checkbox"/>	1.5/2.0	81	12	5.2 x 5.1 x 2	13.3 x 13 x 5.1
120376	Blue	<input checked="" type="checkbox"/>					
120387	White	<input type="checkbox"/>	5	25	22.6	5.2 x 5.1 x 2.9	13.3 x 13 x 7.5
120388	Blue	<input checked="" type="checkbox"/>					
120428	White	<input type="checkbox"/>	Vial	100	17	7.4 x 7.8 x 1.6	18.8 x 19.8 x 4.1




Corrugated Polypropylene 15 mL and 50 mL



Pack of 10

Item No.	Color		Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120391	White	<input type="checkbox"/>	15	36	19.8	5.8 x 5.7 x 4.7	14.7 x 14.5 x 12
120392	Blue	<input checked="" type="checkbox"/>					
120366	White	<input type="checkbox"/>	50	16	31.6	5.8 x 5.7 x 4.7	14.7 x 14.5 x 12
120367	Blue	<input checked="" type="checkbox"/>					



Thin Polypropylene film 1.5–2.0 mL



Pack of 10

Item No.	Color		Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120370	White	<input type="checkbox"/>	1.5/2.0	81	12.1	5 x 4.9 x 1.9	12.5 x 12.4 x 4.9
120371	Blue	<input checked="" type="checkbox"/>					
120372	Lt. Green	<input checked="" type="checkbox"/>					
120373	Purple	<input checked="" type="checkbox"/>					