


Product Highlights


- Chemically-resistant HDPE Coated Wire Racks
- Racks can be used between temperatures of -40°C and 121°C
- Holds all brands and types of conical, round or flat bottom, glass or plastic, tubes and/or bottles
- Square wells and three deck levels help keep tubes secure
- Autoclavable



*HS120089
6x12, 72 Wells

Fits 10-13mm Tubes

Item No.	Color	Tube Array	Wells	Well Diam. mm	L x W x H in	L x W x H cm	UOM
HS120089*	Blue	6x12	72	13	8.1 x 4.1 x 2.6	20.5 x 10.5 x 6.5	1 ea, 48 ea/cs
120759	NEW Blue	9x12	108	13	8 x 6.4 x 2.6	20.4 x 16.1 x 6.5	1 ea, 36 ea/cs



*120761
4x12, 48 Wells

Fits 13-16mm Tubes

Item No.	Color	Tube Array	Wells	Well Diam. mm	L x W x H in	L x W x H cm	UOM
120760	NEW Blue	4x10	40	16	7.9 x 3.3 x 2.4	20 x 8.3 x 5.9	1 ea, 56 ea/cs
120761*	NEW Blue	4x12	48	16	9.4 x 3.3 x 2.4	23.7 x 12.7 x 6.5	1 ea, 56 ea/cs
HS23072	Blue	6x12	72	16	9.5 x 5 x 2.6	24.1 x 18.8 x 6.5	1 ea, 24 ea/cs
120763	NEW Blue	9x12	108	16	9.5 x 7.4 x 2.6	24.1 x 18.8 x 6.5	1 ea, 24 ea/cs

Fits 16-18mm Tubes – 15mL


Item No.	Color	Tube Array	Wells	Well Diam. mm	L x W x H in	L x W x H cm	UOM
HS23215*	Blue	3x5	15	18	5 x 3 x 3	12.7 x 7.6 x 7.6	1 ea, 96 ea/cs



*HS23215
3x5,
15 Wells

Fits 16-20mm Tubes – 15mL

Item No.	Color	Tube Array	Wells	Well Diam. mm	L x W x H in	L x W x H cm	UOM
120764	NEW Blue	4x10	40	20	9.5 x 3.9 x 3.1	24 x 9 x 8	1 ea, 48 ea/cs
120765	NEW Blue	8x10	80	20	9.5 x 7.6 x 3.1	24 x 19 x 8	1 ea, 24 ea/cs



*120768
8x10, 80 Wells

Fits 20-24mm Tubes

Item No.	Color	Tube Array	Wells	Well Diam. mm	L x W x H in	L x W x H cm	UOM
120766	NEW Blue	4x10	40	24	10.2 x 4.2 x 3.2	26 x 10.7 x 8	1 ea, 24 ea/cs

Fits 20-25mm Tubes

Item No.	Color	Tube Array	Wells	Well Diam. mm	L x W x H in	L x W x H cm	UOM
120769	NEW Blue	3x8	24	25	10.7 x 4.1 x 3.2	27.3 x 10.5 x 8	1 ea, 24 ea/cs
120767	NEW Blue	4x10	40	25	11.4 x 4.7 x 3.1	29 x 11.8 x 8	1 ea, 24 ea/cs
120768*	NEW Blue	8x10	80	25	11.4 x 9.1 x 3.1	29 x 23 x 8	1 ea, 12 ea/cs




*HS232461
2x4, 8 Wells

Fits 25-30mm Tubes - 50mL


Item No.	Color	Tube Array	Wells	Well Diam. mm	L x W x H in	L x W x H cm	UOM
HS232461*	Blue	■ 2x4	8	30	5.8 x 3.2 x 3.5	14.7 x 8 x 8.9	1 ea, 96 ea/cs
120770	NEW Blue	■ 6x8	48	30	10.9 x 8.3 x 3.9	27.7 x 20.9 x 9.8	1 ea, 8 ea/cs

Fits 30-40mm Tubes

Item No.	Color	Tube Array	Wells	Well Diam. mm	L x W x H in	L x W x H cm	UOM
120771*	NEW Blue	■ 3x8	24	40	14.1 x 5.5 x 3.9	35.7 x 13.9 x 9.8	1 ea, 12 ea/cs
120772	NEW Blue	■ 6x8	48	40	14 x 10.6 x 3.8	35 x 26 x 9.8	1 ea, 8 ea/cs



*120771
3x8, 24 Wells




*HS120086
5x10, 50 Wells

Fits 23-25mm Bottles - Bijou Bottles

Item No.	Color	Tube Array	Wells	Well Diam. mm	L x W x H in	L x W x H cm	UOM
HS120085	Blue	■ 2x6	12	25	6.4 x 2.5 x 1.3	16.4 x 6.4 x 3.2	1 ea, 96 ea/cs
HS120086*	Blue	■ 5x10	50	25	10.7 x 5.7 x 1.3	27.1 x 14.5 x 3.2	1 ea, 48 ea/cs

Fits 30-33mm Bottles - Universal & McCarthy Bottles

Item No.	Color	Tube Array	Wells	Well Diam. mm	L x W x H in	L x W x H cm	UOM
HS120087*	Blue	■ 2x6	12	33	8.3 x 3.1 x 2.4	21.2 x 8 x 6	1 ea, 72 ea/cs
HS120088	Blue	■ 5x10	50	33	14 x 7.4 x 2.4	35.5 x 18.9 x 6	1 ea, 24 ea/cs



*HS120087
2x6, 12 Wells



*HS120091
4x6, 24 Wells

Fits 55-58mm Bottles - Urine Containers

Item No.	Color	Tube Array	Wells	Well Diam. mm	L x W x H in	L x W x H cm	UOM
HS120091*	Blue	■ 4x6	24	58	15.5 x 10.3 x 2	39.3 x 26.2 x 5	1 ea, 12 ea/cs










Fits 55-60mm Bottles - Centrifuge Bottles





Item No.	Color	Tube Array	Wells	Well Diam. mm	L x W x H in	L x W x H cm	UOM
HS23224*	Blue	■ 2x3	6	60	8.1 x 5.3 x 4.1	20.5 x 13.5 x 10.4	1 ea, 24 ea/cs










*HS23224
2x3, 6 Wells

TUBE STYLE SELECTOR

Example Tube Style									
Example Tube/Bottle/Vial Type	Blood Tubes Culture Tubes Test Tubes	Blood Tubes Culture Tubes Test Tubes	Test Tubes Culture Tubes 15mL Tube	Test Tubes Culture Tubes 15mL Tube	Test Tubes Culture Tubes	Test Tubes Culture Tubes	Test Tubes Culture Tubes 50mL Tube	Test Tubes Culture Tubes	
Well Tolerance	10-13 mm	13-16mm	16-18mm	16-20mm	20-24mm	20-25mm	25-30mm	30-40mm	
Rack Codes	HS120089 120759	120761 HS23072 120763 120760	HS23215	120764 120765	120766	120767 120768 120769	120770 HS232461	120771 120772	

Example Tube Style				
Example Tube/Bottle/Vial Type	Bijou Bottle Scintillation Vial	McCarthy Bottle Universal Bottle Scintillation Vial Culture Bottle Serum Vial	Urine Containers Centrifuge Bottles Culture Bottles Storage Bottles	Centrifuge Bottles Culture Bottles Storage Bottles
Well Tolerance	23-25mm	30-33mm	55-58mm	55-60mm
Rack Codes	HS120085 HS120086	HS120087 HS120088	HS120091	HS23224

TUBE ARRAY EXAMPLES

2x3	2x4	2x6	3x5	3x8	4x6	4x10
						
4x12	5x10	6x8	6x12	8x10	9x12	
