

## Lightweight and easy to use

The RF1000 is a manual pipette controller uniquely designed for one-handed operation. Simply squeeze the main bulb and push the toggle forward to aspirate or back to dispense. Push the blowout bulb to expel any liquid remaining in the pipette.

- Tapered silicone rubber pipette adapter allows a large range of pipette sizes (1-100 mL) Replaceable when worn.
- Accepts common and readily available replacement filters unlike other pipette aids whose filters are custom. This will help reduce total cost of ownership over time.
- Includes one 0.45 µm hydrophobic filter.

Item No.	Description
<b>HS1000RFU</b>	RF1000 Manual Pipette Controller
<b>Replacement parts</b>	
<b>HS5000RFIL</b>	Replacement Filters, pack of 5
<b>120604</b>	Nose Cone with Adapter, Blue

Specifications	
<b>Power:</b>	Manual
<b>Filter:</b>	PTFE/25 mm dia / 0.45 µm filter in PP housing
<b>Pipettes:</b>	1 – 100 mL, glass and plastic
<b>Pipette Cone:</b>	Polypropylene, autoclavable
<b>Dimensions:</b>	180 x 45 x 120 mm



Nose cone assembly only



Filter not included  
in warranty

