

Most versatile rack on the market

- Channels designed for more useful configurations
- Large and Small pegs for all glassware varieties
- One of a kind Drying Basket to collect smaller items

Designed to protect against breakage

- Peg angle allows for better hold - especially large, heavy pieces
- Patent Pending SureGrip Hook securely grips glassware into place

U.S. Patent No. D923,275



Drying Basket

Collects small items and slotted to hold pipets safe while drying



SureGrip Hook

Patent-pending design securely grips glassware into place




Use on counter or wall

- Stands alone on integrated hinged stand
- Mounts easily on wall with included guide and hardware

Easy to clean, durable

Can hold up to 77lbs (34.9 kg)



Lab Drying Rack

Item No.	L x W x H in	L x W x H cm
120791	17.4 x 26 x 5.5	44.1 x 66 x 13.9

Replacement Parts

Item No.	Description
120793	Peg 10 Piece Kit, (5) Small, (5) Large
120792	Drying Basket Kit, 2 ea
120813	Replacement Draining Set, Tube and Elbow

Included:

- 72 Pegs ~ 36 Small (Green) ~ 36 Large (Blue)
- 2 Drying Baskets
- Drainage Tubing
- Mounting Guide & Hardware