

## Deep vessel with sturdy handle designed for convenient handling and pouring of liquids

- Flat bottom ensures stability and prevents tipping
- Notched pour spout design controls liquid pour
- Operating temperature range 0 to 40°C (104°F)
- Made from FDA food grade polystyrene in an ISO Class 8 clean room
- Material biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols and bases
- Recyclable, marked
- Available in 3 Sizes
  - 50mL - 100mL - 200 mL



**Built-in Pour Spout**  
for accurate pouring

**Large Sturdy Handle**  
for ease of movement  
and pouring

200 mL



100 mL



50 mL

**Flat Bottom**  
perfect for placing on  
scales or tabletops  
without tipping



Item No.	Color	Size	Volume mL	L x W x H in	L x W x H cm
120890	White	<input type="checkbox"/> Small	50	3.6 x 1.6 x 1.7	9.2 x 4.0 x 4.3
120891	White	<input type="checkbox"/> Medium	100	4.2 x 3.5 x 1.9	10.7 x 9 x 4.8
120892	White	<input type="checkbox"/> Large	200	5.8 x 4.6 x 2	14.8 x 11.7 x 5.1