

Slide, set, work & store PCR samples

Multi-purpose rack converts from compact workstation to lidded storage solution accommodating individual and strip 0.2 mL tubes

Versatile & Robust Solution

- Polycarbonate lid serves as cover and elevated rack base - just slide open lid pivot lid under rack for compact workstation
- Lid can be used in right or left-handed orientation and is detachable
- Stackable saving valuable bench-top and freezer space
- Multiple colors offered for coding options



US Patent Pending, Community Design No: 004101756-0001-0002

Pack of 5


Item No.	Color	Item No.	Color
120538	Assorted* 	120540	Green 
120539	Blue 	120541	Purple 

*Assorted colors (Purple, Blue, Green, Orange, Yellow)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
0.2	96	5.2	4.9 x 3.6 x 1.5	12.5 x 9.2 x 3.7



Slide, set, work & store



Innovative design

