

Magnetic lab storage series

Easily organize common lab products off your bench top with the MagLab storage system

Vibrant & Functional


- Made of robust ABS plastic in bright translucent colors
- Secure direct to metal with built-in powerful magnets or non-metal surfaces with included self-adhesive mounting plate
- Available individually or as a pack of four

MagRack holds four 50 mL tubes and nine 15 mL tubes





MagWipe holds a standard (110 x 210 mm) Kimwipes® box

MagPette holds two pipettors (fits most popular styles)

MagBox holds miscellaneous accessories



Includes self-adhesive metal plates that adhere to any surface

Item No.	Description	Color	L x W x H in	L x W x H cm
HS23500	MagLab Series, Pack of 4	Assorted 	—	—
HS23502	MagWipe	Green 	5.1 x 3.5 x 3.8	13 x 9.1 x 9.8
HS23503	MagPette	Purple 	5.1 x 4.4 x 2.9	13.1 x 11.2 x 7.6
HS23504	MagBox	Orange 	5.1 x 3 x 4	13 x 7.8 x 10.2

*Assorted colors: Teal, Green, Purple and Orange