

Microtube Racks

Durable, Simple to Use

- Made from reinforced nylon offering superior chemical and heat-resistance over most other plastic based racks
- Handles provide easy and safe transportation
- Stable design can be used in water baths
- Alphanumeric identification for quick location and retrieval of tubes



Item No.	Color	Tubes mL	Wells	Well Dia. mm	Packaging	L x W x H in	L x W x H cm
120241	Blue	■ 1.5/2.0	24	11.1	8	5.6 x 3.5 x 2	14.3 x 9 x 5.2
120243			96		4	8.3 x 5.2 x 2.1	21 x 13.1 x 5.4

Microtube Racks 80-well

Popular fraction-collector-plate style rack

Compact 5x16 array makes filling and transfer easy

- Frosted labeling area on bottom of rack
- Embossed index numbers for quick sample ID and retrieval
- Made of polypropylene



Pack of 5

Item No.	Color	Part No.	Color
HS29025A	Natural □	HS29025E	Yellow ■
HS29025B	Blue ■	HS29025F	Orange ■
HS29025C	Green ■	HS29025G	Assorted* ■ ■ ■ ■ ■
HS29025D	Pink ■	HS29025J	Red ■

*Assorted colors (Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
1.5/2.0	80	11	8.9 x 2.6 x 1	22.5 x 6.7 x 2.8





Microtube storage boxes

Safe microtube storage

Low profile boxes are available in two sizes to accommodate standard 0.5 mL and 0.2 mL microtubes and vessels of a similar size

- Partitions included in the boxes
- Will fit into most standard freezer racking



Pack of 10

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120220	White	□ 0.2	169	9.5	5.2 x 5.2 x 1.1	13.3 x 13.3 x 2.9
120221	White	□ 0.5	100	12	5.2 x 5.2 x 1.6	13.3 x 13.3 x 4.1

Cryogenic Vial Boxes & Partitions

For storage of small samples and tubes

- Made from cardboard
- Colors allow for easy identification and sorting purposes
- All box footprints are 13.3 x 13.3 cm, select from two different heights



Vent Holes

allow faster cooling and removal of LN₂

Cardboard Cryogenic Vial Boxes

Pack of 12

Item No.	Color	Tubes mL	L x W x H in	L x W x H cm
HS2860A	White	□	5.2 x 5.2 x 2	13.3 x 13.3 x 5
HS2860B	White	□	5.2 x 5.2 x 3	13.3 x 13.3 x 7.5
HS2860C	Assorted	■ ■ ■		
HS2860CB	Blue	■		
HS2860CR	Red	■	5.2 x 5.2 x 2	13.3 x 13.3 x 5
HS2860CY	Yellow	■		

Partitions

Partitions must be ordered separately

Pack of 12

Item No.	Array	Tubes mL	Wells	Well Dia. mm
HS2860E	10 x 10	1.5/2.0	100	11.5
HS2860H	9 x 9	1.5/2.0	81	12
HS2860F	8 x 8	1.5/2.0	64	15
HS2860G	7 x 7	1.5/2.0	49	17.5