

## Mini Vortexer

### Small footprint with maximum power

The Mini Vortexer by Heathrow Scientific offers immediate vortex shaking for a variety of tube sizes. A simple touch start operation enables exceptional mixing performance with one hand


**5**  
Year  
Warranty

#### Robust Design

- Low profile design and small footprint saves valuable space
- Made of ABS and silicone materials that resist most commonly used chemicals
- Stable platform and elastomer feet ensure unit stays stable during operation

#### Powerful Function

- Rapid vortex action with low activation force
- Quiet operation and low vibration



**Non-slip Feet**  
ensures the unit stays  
in place while mixing

**5**  
Year  
Warranty

**Did you know?**

Heathrow Scientific has strict standards that include 3<sup>rd</sup> party plant reviews and 100% multi-level product inspections. These actions result in a 99.7% reliability on our equipment and 99.99% reliability on our laboratory supplies.

**Flexible**

Can be used with tubes ranging from 2.0 mL to 50 mL



U.S. and E.U. Utility Patent Pending  
U.S. Design Patent No. D814,045

Item No.	Color	L x W x H in	L x W x H cm
120567	Grey/Purple	5.1 x 4.3 x 3	13 x 11 x 7.5
120598	Grey/Blue		
120827	Black/Black		

**Specifications**

<b>Speed</b>	4,500 rpm
<b>Orbit</b>	6 mm
<b>Capacity</b>	Up to 50 mL tube



**Replacement parts**

120617	Mini Vortexer Silicone Cover, set of 2
120814	Mini Vortexer Universal Power Adapter

**Package Includes**

- Mini Vortexer
- 2 Silicone Caps
- 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs

