




Sliding Lock
provides secure storage with easily identified lock and unlock symbols.

Slide Rail Design
ensures a secure fit allowing easy placement and removal of slides

Centrally Positioned Numbering
enables simple slide organization with easy viewing of location numbers



100-Place Premium Plus Slide Box

Modern design for maximum functionality

Optimal Organization

- Organize by "project" and "date" molded in on lid
- Inventory card index corresponds to slide numbers and is removable for trouble-free writing on flat surface
- Vibrant translucent colors are easy to identify and enable coding options

Functional Design

- Flexible storage, store vertically or horizontally to efficiently use space
- Molded feet allow secure horizontal stacking of multiple units
- Unique side storage can accommodate a pen, desiccant or other small items needed for working
- Hinged lid smoothly opens 180° sitting flat on surfaces and prevents misplacement

Safe & Secure

- EVA foam lines the bottom of the box helping to hold slides securely in place and protect against damage
- Stylish wave surface texture is easy to hold for safe sample transport

100-Place Premium Plus Slide Box

Pack of 5

Item No.	Color	Item No.	Color
120576	Blue	120579	Red
120577	Orange	120592	Purple
120578	Green	120593	Assorted*

*Assorted colors (Purple, Blue, Green, Red, Orange)

Lining	Lock	L x W x H in	L x W x H cm
EVA Foam	Sliding Latch	8.2 x 7.7 x 1.3	20.7 x 19.5 x 3.2

Replacement Inventory Card

Pack of 10

Item No.	L x H in	L x H cm
120618	8.5 x 11	21.6 x 28

Premium Microscope Slide Boxes

Slide boxes are perfect for storage or transportation of standard size slides (76 x 26 mm and 3 x 1 in)

Organize & Identify

- Lining and slots are numbered to correspond to the slide inventory sheet
- Organize by "Project" and "Date" with pre-molded labeling area
- Available in multiple colors for project coding

Durable & Safe

- ABS plastic boxes available with either cork or foam lining
- Designed to securely store slides in grooved slots
- Molded feet secure stacking for easy transport and storage



25-Place Slide Box

Easily access slides with friction-fit cover



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15989A	Blue <input checked="" type="checkbox"/>	Cork	Friction Fit	5.6 x 3.5 x 1.4	14.1 x 8.8 x 3.5
HS15989B	Green <input checked="" type="checkbox"/>				
HS15989C	Red <input checked="" type="checkbox"/>				
HS15989E	White <input type="checkbox"/>	Foam	Friction Fit	5.6 x 3.5 x 1.4	14.1 x 8.8 x 3.5
HS15989M	Blue <input checked="" type="checkbox"/>				
HS15989P	Red <input checked="" type="checkbox"/>				
HS15989R	White <input type="checkbox"/>				

50-Place Slide Box

Thumb-latch lock mechanism



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15996A	Blue <input checked="" type="checkbox"/>	Cork	Thumb-Latch	8.2 x 3.4 x 1.4	20.8 x 8.6 x 3.5
HS15996B	Red <input checked="" type="checkbox"/>				
HS15996C	White <input type="checkbox"/>				

100-Place Slide Box

Closes securely with a rust-resistant nickel-plated clasp and pin hinges



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15994A	Blue <input checked="" type="checkbox"/>	Cork	Nickel-plated Clasp	8.2 x 6.9 x 1.3	20.8 x 17.5 x 3.4
HS15994B	Green <input checked="" type="checkbox"/>				
HS15994C	Red <input checked="" type="checkbox"/>				
HS15994D	Yellow <input checked="" type="checkbox"/>				
HS15994E	White <input type="checkbox"/>				
HS15994F	Grey <input type="checkbox"/>				
HS15994G	Black <input checked="" type="checkbox"/>	Foam	Nickel-plated Clasp	8.2 x 6.9 x 1.3	20.8 x 17.5 x 3.4
HS15994M	Blue <input checked="" type="checkbox"/>				
HS15994P	Red <input checked="" type="checkbox"/>				
HS15994R	White <input type="checkbox"/>				

100-Place Tall Slide Box

Additional height for 3 x 2 in slides
Stackable with 100-Place Premium Slide Boxes



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15979A	Blue <input checked="" type="checkbox"/>	Cork	Nickel-plated Clasp	8.2 x 6.9 x 2.4	20.8 x 17.5 x 6