

Droplet™ Corrugated Sample Storage Tray

Corrugated storage trays can be used with the included dividers for specimen containers or remove them for general lab product storage and transport

- Made from corrugated polypropylene
- Moisture resistant, significantly reducing the chance of mold
- Stackable with other trays
- Pack includes: 10 trays and 10 dividers


Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120375	White <input type="checkbox"/>	24	58	16 x 10.2 x 3.5	40.6 x 26 x 8.9



Use for storage or transport of laboratory products



Removeable inserts



Trays are supplied "flat-packed"
Reducing the amount of space needed for storage



Lightweight and Sturdy
to hold up to 24 full specimen containers

Droplet™ Sample Storage Tray

Sample organization for the busy lab

Hard wearing storage trays that can be used for specimen containers with the insert provided, or for general lab products without the insert

- Lightweight and sturdy to hold up to 24 full specimen containers (Specimen containers not included)
- Stackable with other trays to support as many as four high
- Can be easily cleaned with common disinfectants and germicides
- Bin and tray are made from polystyrene



Removable Inserts
allow outer shell to be used for storage or transport of laboratory products

Integrated Label Holder




Use for storage or transport of laboratory products

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120172	Blue/White	24	58	18.9 x 13 x 3.3	48 x 33.1 x 8.4

Coated Wire Rack - Urine Containers

Durable racks for transport or sample handling

- Chemically-resistant, HDPE-coated steel wire rack
- Fits bottles with 55-58 mm diameters
- Mesh bottom support



Item No.	Color	Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120091	Blue	6 x 4	24	58	15.5 x 10.3 x 2	39.3 x 26.2 x 5