

## Pipettes never looked so good on your benchtop

Stylish pipette stand holds and stores up to 8-single channel, 4-multi channel, or a combination of both pipettes

### Versatile



- Patent pending design enables this stand to securely hold virtually every known brand of pipette in the marketplace
- Modular upper holder inserts accommodate the fit of different pipette widths
- Assorted color inserts allow for color coding configurations

### Durable Utility

- Stable weighted base allows for easy one-handed removal and replacement of pipettes
- Rotating head makes pipette selection fast and easy
- Moving the stand while keeping pipettes in place is simple with a quick grip handle allowing a user to grab and go



Available in solid or a combination of colors



Assorted color option allows for color coding pipette location



Item No.	Upper Holder Color	Quantity	Lower Holder Color	Quantity
120480	Blue/Green	8 ea/color	Blue	8 ea
120481	Purple/Clear	8 ea/color	Purple	8 ea
120499	Assorted/Clear*	1 ea/color 8 clear	Clear	8 ea

\*Assorted colors (Purple, Blue, Green, Grey, Red, Orange, Yellow and Clear)

Base Width in	Top W x H in	Base Width cm	Top W x H cm
6.2	9 x 16.1	15.8	22.9 x 40.8