

A truly green solution to sample weighing

Made from FSC® Certified (FSC® Mix Credit) paper (Eco-IN) and when disposed of in the standard waste stream they will naturally biodegrade (Eco-OUT)

- A cost-effective solution to reduce plastics in the lab
- Made of grease-proof paper that meets FDA 21 CFR 176.170 (components of paper in contact with aqueous and fatty foods) and with 21 CFR 176.180 (contact with dry foods)
- Control contents like a polystyrene weigh boat because they are easier to pour
- An alternative to regular weighing papers because they are easier to handle

Item No.	Color	Size	Packaging	L x W x H in	L x W x H cm
120312	Brown <input checked="" type="checkbox"/>	Small	Pack of 500	1.4 x 1.4 x 0.6	3.5 x 3.5 x 16
120314	White <input type="checkbox"/>	Small	Pack of 500	1.4 x 1.4 x 0.6	3.5 x 3.5 x 16
120313	Brown <input checked="" type="checkbox"/>	Medium	Pack of 500	1.9 x 1.9 x 0.6	4.8 x 4.8 x 14
120315	White <input type="checkbox"/>	Medium	Pack of 500	1.9 x 1.9 x 0.6	4.8 x 4.8 x 14
120626	Brown <input checked="" type="checkbox"/>	Large	Pack of 250	3.8 x 3.8 x 1	9.5 x 9.5 x 2.5
120195	White <input type="checkbox"/>	Large	Pack of 250	3.8 x 3.8 x 1	9.5 x 9.5 x 2.5

smartBoats™ is a trademark of LevGo Inc.



Functional in the Lab

Available in a small and medium in both brown for light colored samples and white for dark samples

